

FAM109B-family with sequence similarity 109, member B Gene

Reference: PTXBC104175

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	IPIP27B; Ses2; sesquipedalian-2; 27 kDa inositol polyphosphate phosphatase interacting protein B; family with sequence similarity 109 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 109, member B
GENE ID	150368
ORIGIN SPECIES	Human
ACCESSIONS	BC104175
VECTOR	pDONR223
NCBI SYMBOL	FAM109B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgaacgagaggagttagccactatgcactcagcactccccagcggaccacatggcttctcgcgacctggggggccaggacccccaccgacccccagtgccactggccgaagatctggtttgtctcaagggaacctgctattctctttgagagtcgcgagg
gccccggcccccactgagcctgtgtgtgctggaaggtgcacagtggaactggccgaggtccccgtgccccgaggagtttgcctttgccatctgctttgatgccccctggagtgccccacactgtgctgcccgcagaaggccggcgccccaggagcctgggtgaaggtgctgtccgg
gcaagctttggctacatgcgcctgtgtgtacgcgagttggagagccagttgcaggacgcacgcccagagcctggctttgcaacgccgtcatctctggaagtctgttgccagccgctgtaagccccaggctcctaaccaccgagctgcgggctggagaalggccactgcctctccaagga
cagcagccctgtgggctgtgtgaagaagcgggcagcaggtctgcagggtgggggtgtgctgagtgaggagctgcagggccctgccagctctccttagcgaaggggcagagccctgtgtccctgagacctctgtctctaccctgcatgactggtatggccaggagatctggag
ctgcggcagtggtggcagaagaggcccccaggggagccactcaaaatgtgaggaacagatagggcctaa