

LOC285194-hypothetical LOC285194 Gene

Reference: PTXBC104179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LINC00902; LSAMP-AS1; LSAMP-AS3; NCRNA00295; tumor suppressor candidate 7 (non-protein coding)		

For research use only.

GENE DESCRIPTION	hypothetical LOC285194
GENE ID	285194
ORIGIN SPECIES	Human
ACCESSIONS	BC104179
VECTOR	pDONR223
NCBI SYMBOL	LOC285194
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaattcagagaacatgacttcctattatccttctcaaagttatgctgcagaggaaagaagcatacatcttttaccaccaggaaaccccaagcatctattaccataatagccatgggaaacagaaggcacctcaataaagggtgggaaaaaatgaaagaaatggcttggcctgtg
cctgtttgaccttgagagatactttttgcaagaattgtaagttttggcccaaaaagtgtcggtcttttatccttctctgtggagccaaactgcaaaccaggatga