

LL22NC03-75B3.6-KIAA1644 protein Gene

Reference: PTXBC104183

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	KIAA1644 protein
GENE ID	85352
ORIGIN SPECIES	Human
ACCESSIONS	BC104183
VECTOR	pDONR223
NCBI SYMBOL	LL22NC03-75B3.6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccagttgtggccagcagtccttgaacgtgctcgccgtctctctcattgtgttttctgcagtccttctgtcacatttcgggtctgtgaaccatacacagaccacaaaggccgctaccatttggcttccactgccccggctctcggaacaagacattcatctgtgttcaccataa
caacacgggtcttcaataactgctgcaacgagagggagttccaggcgggtgatgcaggcgaacctcacggccagctccgaggggttacatgcacaacaattacaccgccctgttgggagtggtgatctatggattttcgtgttgatgctgctggttctggacctttgtattactcggcaatgaact
acgacatctgcaaggctctacctggcaggtgggcatccaaggacgatggatgaacaggacccccggcgtgggggaaccccgctcgggccccctcggccgggtcagcgggccccacagccgcagcctccccaggcccgctgccacaagccccacaggccgtgcacacattg
cggggagatgtcacagcccaccgtgatgacctccagagttcgtctgcctga