

RIT1-Ras-like without CAAX 1 Gene

Reference: PTXBC104186

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GTP-binding protein Rit1; NS8; RIBB; RIT; ROC1; GTP-binding protein Roc1; Ric-like, expressed in many tissues; ras-like protein expressed in many tissues; ras-like without CAAX protein 1; Ras like without CAAX 1		

For research use only.

GENE DESCRIPTION	Ras-like without CAAX 1
GENE ID	6016
ORIGIN SPECIES	Human
ACCESSIONS	BC104186
VECTOR	pDONR223
NCBI SYMBOL	RIT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattctggaactcggccagttggttagctgctgtagcagccccgctgggctctcagggagtagcaaaactagtgtgctgggtgctggtggtgtaggggaagagtgccatgacctgcagttcatcagccaccgattccagaagatcatgcccaccattgaagatgcttataagatcag
gatccgtattgatgatgagcctgccaatctggacattttggatacagctggacagggcagagtttacagccatgcgggaccagtatatagggcaggagaagggtttatcatctgttactctatcacggatcgtcgaagtttccatgaagttcgtgagtttaaacagcttattatcgagtcggacgt
actgacgatacaccctgtggtctgtggtggaacaagtcagacctcaaacagctaaagacaggtcaccaagggaagaggtggccttggcccgagaattcagctgtccctttttgagacatctgctgcataccgctactatatgtgatgttttccatgccctgtacgggagatacgttaggaaa
gaaaaggaggcagtagtggccatggagaaaaaatcaagccaaaaacagtgtagtgaagaggctaaatcaccattccggaagaagaagattcagtaacttga