

NEDD8-neural precursor cell expressed, developmentally down-regulated 8 Gene

Reference: PTXBC104201

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ubiquitin-like protein Nedd8; NEDD-8; neddylin; neural precursor cell expressed, developmentally down-regulated 8		

For research use only.

GENE DESCRIPTION	neural precursor cell expressed, developmentally down-regulated 8
GENE ID	4738
ORIGIN SPECIES	Human
ACCESSIONS	BC104201
VECTOR	pDONR223
NCBI SYMBOL	NEDD8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctaattaaagtgaagacgtgaccggaaggagattgagattgacattgaacctacagacaagggtggagcgaatcaaggagcgtgtggaggagaaaagggaatccccccacacagcagaggctcatctacagtggcaagcagatgaatgatgagaagacagcagctgattac
aagattttagggtgttcagtccttcacctggtgtgtgctctgagaggagggtgtgtcttaggcagtga