

DPM3-dolichyl-phosphate mannosyltransferase polypeptide 3 Gene

Reference: PTXBC104202

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDG10; dolichol-phosphate mannosyltransferase subunit 3; DPM synthase complex subunit 3; DPM synthase subunit 3; MPD synthase subunit 3; dolichol-phosphate mannose synthase subunit 3; dolichyl-phosphate beta-D-mannosyltransferase subunit 3; dolichyl-phosph		

For research use only.

GENE DESCRIPTION	dolichyl-phosphate mannosyltransferase polypeptide 3
GENE ID	54344
ORIGIN SPECIES	Human
ACCESSIONS	BC104202
VECTOR	pDONR223
NCBI SYMBOL	DPM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctcgtggcgccgcttcggttgagtttggtccgctttctctctgctcctcagggagcattgctctctctcgcagtgaccatgacgaattagcgcagtggtcttgggactagcgcgtctgggtccacctgggtggccctgaccacgggagccttgggcctggagctgccctt
gtctgccagggaagtctgtgcccactgcccgcctactgtgtgtgtccggtgtgctatgcccctgggcactgtgggtatctgtgtggccacttttcatgactgcagagacccgcacgcagctgcagagccagatacaggagcccagccgacttagcccgaggggctgcgc
ttctga