

NCRNA00085-non-protein coding RNA 85 Gene

Reference: PTXBC104208

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NCRNA00085; LET7EH; LINC00085; SPACA6P; sperm acrosome membrane-associated protein 6; BACHELOR-like protein; Let-7e host; long intergenic non-protein coding RNA 85; sperm acrosome associated 6, pseudogene; sperm acrosome associated 6		

For research use only.

GENE DESCRIPTION	non-protein coding RNA 85
GENE ID	147650
ORIGIN SPECIES	Human
ACCESSIONS	BC104208
VECTOR	pDONR223
NCBI SYMBOL	NCRNA00085
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctgtggtctgcccagtgccgtccctgtgccctgtggccctggctgtctcaggggtgcccgctgggctgtctcctctgtctcacaactactctgagcgcctccgcatctgccagatgtttgtggatgaggagcccaagcttgaagagtgtgaggagccttcacggccgccttcaggcctctgacaccgaaatcagtgaggagaccatccacacttcatcagtgctctgggaaggtgcagagggagggcaggagagggccagaggtcaggtgagggacagacagagagaacagtcagaggagaaggtctaaagaccatgagaacaacagagacttagggacagagacacagacaggggaagacagcagggcaaaagactcagagaggggaggtggagagtcagagaggggaagatggagactcagagaggggaggtggagactcagagagaggaagatggagactcggaaagatggagactcaggagatggagagtcagagagggaggtggacactcggaggtggagagtcaggaggtggagactatagaaggggaggtggagagtcagagaggtggagactggagaggggaatag