

ENDOGL1-endonuclease G-like 1 Gene

Reference: PTXBC104212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ENDOGL1; ENDOGL2; ENGL; ENGL-a; ENGL-b; ENGLA; ENGLB; nuclease EXOG, mitochondrial; endo G-like 1; endo/exonuclease (5'-3'), endonuclease G-like; endonuclease G-like 1; endonuclease G-like 2; exo/endonuclease G		

For research use only.

GENE DESCRIPTION	endonuclease G-like 1
GENE ID	9941
ORIGIN SPECIES	Human
ACCESSIONS	BC104212
VECTOR	pDONR223
NCBI SYMBOL	ENDOGL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccagcaggaatacaaatctcaagtaaacccatggctgaacctttacatttgaacattgtgcctcaggatttgaataaattctgataattggaacagaatagaatgtactgtcgagagctgacagaaaggtttgaagatgtttgggtgtatctgggcctttgacctacacag
actagaggcgtggaagaaaatagtttagttaccaggtgattggcgaggacaacgtggcagtcaccttataaggttaacctggccgcagaaagctcagatctacgaaccactggcgctaggggcctttgtgtaccatgaagccatcggttcagccccagtttaactga
attccaagtgcagcctcaggacacagagattgtcaggactgggtgtttttctcatttggatagaactagtgataccggaatatctgctctgtggacacctgtaagctcctggaattccaggagttcacctgtacttgagtacaagaagattgaaggagcccgatcagtgctcagactgaa
aagatcatgaaaactgaagaatgcagagattgaaccagatgattactttatgagtcgctatgagaagaagctagaagaactcaagctaaggagcagtcaggaaccagataagaagccatcctag