

FLJ25328-hypothetical LOC148231 Gene

Reference: PTXBC104221

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC148231
GENE ID	148231
ORIGIN SPECIES	Human
ACCESSIONS	BC104221
VECTOR	pDONR223
NCBI SYMBOL	FLJ25328
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatagttggccatggctctttgatagcaattttgtctgcccgtagaaattatccttgaaccagaaaaaaatcccaattcaacactccattttgtgaagagctggtccgaattctgggtcacttttggaaaaatatctgttggtaaacagtgcctcatctgtgatcgaagtctctactg
gaattggatagagatggtgaatttcagccagagcgcaccaagaactcgtgtcccgccgtcctgaacgccatcagccgcgcacctccatctccctcccccagagtcctctgctcccccgtgactctggcctcctgcagtggttctgtggccacgtcaggaagtgtggagacgtggca
ttgcaccaaacaccaggcggccagccgtgtggtcctccgcaaatcaggctcttctgtgcatggaggtggagtcctctggaggcagccgtggccctgttctcgggcaggtggagcatggagacgactgggaactcactagacacgcacccagccggttcaagggtgtcgcgc
gtgcaccagcttgcaggagaaagcccaacgcctcgaagcacgaaggcctccctcccgacttga