

CCL4-chemokine (C-C motif) ligand 4 Gene

Reference: PTXBC104226

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACT2; AT744.1; G-26; HC21; LAG-1; LAG1; MIP-1-beta; MIP1B; MIP1B1; SCYA2; SCYA4; C-C motif chemokine 4; G-26 T-lymphocyte-secreted protein; MIP-1-beta(1-69); PAT 744; SIS-gamma; T-cell activation protein 2; chemokine (C-C motif) ligand 4; lymphocyte activ		

For research use only.

GENE DESCRIPTION	chemokine (C-C motif) ligand 4
GENE ID	6351
ORIGIN SPECIES	Human
ACCESSIONS	BC104226
VECTOR	pDONR223
NCBI SYMBOL	CCL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctctgctgactgtctctctcatgtagctgcttctgctcagcgctcagaccaatgggtcagaccccccacgcctgctgctttcttacaccgcgaggaaagcttctcgaactttgtgtagattactatgagaccagcagcctctgtcccagccagctg
tggatattccaaaccaaagaagcaagcaagtctgtgctgacccagtgagacctgggtccaggagtacgtgtatgacctggaactgaactga