

ATP6V1F-ATPase, H⁺ transporting, lysosomal 14kDa, V1 subunit F Gene

Reference: PTXBC104230

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ATP6S14; VATF; Vma7; V-type proton ATPase subunit F; ATPase, H ⁺ transporting, lysosomal 14kDa, V1 subunit F; ATPase, vacuolar, 14 kD; H(+)-transporting two-sector ATPase, 14kD subunit; V-ATPase 14 kDa subunit; V-ATPase F subunit; V-ATPase subunit F; adeno		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal 14kDa, V1 subunit F
GENE ID	9296
ORIGIN SPECIES	Human
ACCESSIONS	BC104230
VECTOR	pDONR223
NCBI SYMBOL	ATP6V1F
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggagggttaagctcatgcagtgatcggagacgaggacacggtgactggttctctgctggcgccataggggagcttaacaagaaccgccatcccaatttctggtgtggagaaggatacaaccatcatgagatcgaagacactttccggcaatttctaaaccgggatga
cattggcatcatcctcatcaaccagttacatgcagagatgggtcgccatgcctggacgcccaccagcagtcctatcccgctgtctctggagatccctccaaggagcaccatagacgccgccaaggactcctctgcagggccagggcatgttctactgccgaagacctgcct
ag