

LCE1B-late cornified envelope 1B Gene

Reference: PTXBC104232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LEP2; SPRL2A; late cornified envelope protein 1B; late envelope protein 2; small proline rich-like (epidermal differentiation complex) 2A; small proline-rich-like epidermal differentiation complex protein 2A; late cornified envelope 1B		

For research use only.

GENE DESCRIPTION	late cornified envelope 1B
GENE ID	353132
ORIGIN SPECIES	Human
ACCESSIONS	BC104232
VECTOR	pDONR223
NCBI SYMBOL	LCE1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcctgccagcagaaccagcagcagtcgagcccccctcccaagtcaccccccaagtcgctcaccctagatgcccccaagtcgcccccaagtcgctcagctctctctgctgagtcagctccggaggtgctgtggtccagctctgggggaagctgtg
gctccagctctgggggtgctgagttctgggggaggtggctgctgagccaccacagcgccgtaggtccactgccacagccccagagctctggtgctgagccagccctccggggctccagctgctgtggaggaggagtgccagcactctggaggtgctgctaa