



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ASIP-agouti signaling protein, nonagouti homolog (mouse) Gene

Reference: PTXBC104238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGSW; AGTI; AGTIL; ASP; SHEP9; agouti-signaling protein; agouti signaling protein, nonagouti homolog; agouti switch protein; nonagouti homolog; agouti signaling protein		

For research use only.

GENE DESCRIPTION	agouti signaling protein, nonagouti homolog (mouse)
GENE ID	434
ORIGIN SPECIES	Human
ACCESSIONS	BC104238
VECTOR	pDONR223
NCBI SYMBOL	ASIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgtcaacccgcttactctgtgccaccctctgctgttctctgcttcttactctgccaacagccacctgccacctgaggagaagctcgcgagatgacaggagcctctgagaagcaactctctgtgaacctactggatgtccctctgtctctatgtggcgctgaacaagaatatcaaacagatcggcagaagaagcagacagaaaagaagatcttctaagaaggagcgttcgatgaagaaagtgtgtcgcgcccgaccgccacctatctgcgcctctgctgtgccaccgcacaacagctgcaagccgccgccaccgcctctgctgcgaccgtgctgcctctgacagtgcgccttcttcgcacgcccctctgcctctgcctctcaactctcga