

C9orf27-chromosome 9 open reading frame 27 Gene

Reference: PTXBC104241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C9orf27; EST-YD1; long intergenic non-protein coding RNA 474		
For research use only.			

GENE DESCRIPTION	chromosome 9 open reading frame 27
GENE ID	58483
ORIGIN SPECIES	Human
ACCESSIONS	BC104241
VECTOR	pDONR223
NCBI SYMBOL	C9orf27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctctgttagtgattcacaagcaacctgtactctcaagaagaagatgataagatgaaggagatcttagaaccagcaactggggctcttcttcagtgaagaaagtgggtgatgcagactcatccatccatgaatgattgcagggatgtgacctatgtaatgaactgctactattgctcatcatctcttcaatga