

DYNLL1-dynein, light chain, LC8-type 1 Gene

Reference: PTXBC104245

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DLC1; DLC8; DNCL1; DNCLC1; LC8; LC8a; PIN; hdlc1; dynein light chain 1, cytoplasmic; 8 kDa dynein light chain; cytoplasmic dynein light polypeptide; dynein, cytoplasmic, light polypeptide 1; protein inhibitor of neuronal nitric oxide synthase; dynein ligh		

For research use only.

GENE DESCRIPTION	dynein, light chain, LC8-type 1
GENE ID	8655
ORIGIN SPECIES	Human
ACCESSIONS	BC104245
VECTOR	pDONR223
NCBI SYMBOL	DYNLL1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcgaccgaaaggccgtgatcaaaaatcgagacatgtcgaagagatgcaacaggactcggtagtgctgctactcaggcgttgagaaatacaacatagagaaggacattcggtctcatcaagaaggaaattgacaagaagtacaatcccacctggcattgcatcgtggga
ggaactcggtagttatgtgacacatgaaaccaaacacttcactctctacctggccaagtggccattctctgttcaaatctggttaa