

LOC150185-hypothetical LOC150185 Gene

Reference: PTXBC104251

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	long intergenic non-protein coding RNA 895		
For research use only.			

GENE DESCRIPTION	hypothetical LOC150185
GENE ID	150185
ORIGIN SPECIES	Human
ACCESSIONS	BC104251
VECTOR	pDONR223
NCBI SYMBOL	LOC150185
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggcccgagatgtgtgccctctccttgctcctgaccgctgactcgttattccaaatagcagccgaaaccatcgctcttcaggtgcagcttttctctgcctctctatgtggaagcagcgcagcgtggactccagccaaactccctcgattggattgcagctctctggtctctgtcctg
tgacctgggcagtttacctcctctatgtctcagtttccctcagctgtgcaaaggggttactgacctgcttaatgggtgggtgtgtgaagagggaatgagttaa