

IL10-interleukin 10 Gene

Reference: PTXBC104252

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSIF; GVHDS; IL10A; TGIF; interleukin-10; T-cell growth inhibitory factor; cytokine synthesis inhibitory factor; interleukin 10		

For research use only.

GENE DESCRIPTION	interleukin 10
GENE ID	3586
ORIGIN SPECIES	Human
ACCESSIONS	BC104252
VECTOR	pDONR223
NCBI SYMBOL	IL10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacagctcagcactctgttgcctggtcctcctgactgggtgagggccagccaggccaggccagcagctctgagaacagctgcacccactccaggccaactgccttaacatgcttcgagatctccgagatgccttcagcagagtgaagactttcttcaaatgaaggatcagctggacaacttgtttaaaggagctctgtgaggactttaagggttacctgggttgccaagcctgtctgagatgatccagttttacctggaggaggtgatgcccaagctgagaaccaagaccagacatcaaggcgcatgtgaactccctgggggagaacctgaagacctcaggctgaggctacggcgctgcatcgatttctccctgtgaaaacaagagcaaggccgtggagcagggtgaagaatgcctttaataagctccaagagaaaggcatctacaagccatgagtgagtttgacatcttcatcaactacatagaagcctacatgacaatgaagatacgaaactga