

## C8orf83-chromosome 8 open reading frame 83 Gene

Reference: PTXBC104255

|                |                                                                                                                                                                        |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | C8orf83; PRO0845; UPF0599; triple QxxK/R motif-containing protein; protein TRIQK; triple repetitive-sequence of QXXK/R protein homolog; triple QxxK/R motif containing |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 8 open reading frame 83                                                                |
| <b>GENE ID</b>          | 286144                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC104255                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C8orf83                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggtaacctcagataaccattttaaatctgctttgttgcctatatactctgtctcctaatgattcagcacaacaacattccagagagaaaaaacaagagttctggaggctgctcacattccttgcttatggcctcctcttccaacgtcaaggccagcaacattgcagcaatctgacttatt  
tgtttgatgtctcctgaccacatctggaagggttatctacttttgaagacacatccagataattcaggacaatctgtccattttgaggttctcattga