

PKIG-protein kinase (cAMP-dependent, catalytic) inhibitor gamma Gene

Reference: PTXBC104257

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PKI-gamma; cAMP-dependent protein kinase inhibitor gamma; protein kinase (cAMP-dependent, catalytic) inhibitor gamma		

For research use only.

GENE DESCRIPTION	protein kinase (cAMP-dependent, catalytic) inhibitor gamma
GENE ID	11142
ORIGIN SPECIES	Human
ACCESSIONS	BC104257
VECTOR	pDONR223
NCBI SYMBOL	PKIG
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggaggtcgagtcctctactcggactcatctcgtgtgacccgacagccgtcggaatgcgtccctgacatccagggagactcagaggctgtgagcgtgagggaagctggctggagacatggcgagctggcactcagggggcagaaggacaggtggagggaagcggcc
cagacaaggaaagctggcaaccagccccagagcagcgtatggaccacctgtcttga