

CXorf27-chromosome X open reading frame 27 Gene

Reference: PTXBC104428

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CXorf27; HIP17; huntingtin-interacting protein M; huntingtin yeast partner M; huntingtin interacting protein M		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 27
GENE ID	25763
ORIGIN SPECIES	Human
ACCESSIONS	BC104428
VECTOR	pDONR223
NCBI SYMBOL	CXorf27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcagagaaaaaactgtgaagaactgtctacaataacaaccagactcaagacccttctagaatgagctacaggtccctaggagcttctgtgaccgcgttgcaagatgaacgagacgtccaaagccagagttcctccacaataataaccctctaccttgctcgattgccttgctg
actacatcatggagcgggtaggcttgaggccagcaacaatggcagtatgcgaacacttcacaagataggagagagaagtggacaacaaccgtgagcccccacagcgtgagagtgtgtgactgcctttttgtttgatgatgcccgaatccaggagaatgactga