

DKFZp566H0824-hypothetical LOC54744 Gene

Reference: PTXBC104430

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC54744
GENE ID	54744
ORIGIN SPECIES	Human
ACCESSIONS	BC104430
VECTOR	pDONR223
NCBI SYMBOL	DKFZp566H0824
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacagcacatagagctggaggggtggccagcagcaaccaggagagaaggcgttcagggtgtctgcaccaaggagatcctgggtcagagccaggacatgtgggaaaaaccattggtccctttaccagggaaggactcaagaagccgaacgtggtaggagatggaggctccc
aagccctggcagacgttacagcagcgtgagcacttcaggaccaggaactccaagggtccagctctggaacgccctctctgaccttgggttcctcatcagcagagggtctctgggtctatgggtcaagctcgaagtcaccttaagagaaatctacaccttctgtctccttcattgt
ga