

HCG27-HLA complex group 27 Gene

Reference: PTXBC104432

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bCX101P6.9; bPG299F13.9; bQB115I13.2; HLA complex group 27 (non-protein coding)		

For research use only.

GENE DESCRIPTION	HLA complex group 27
GENE ID	253018
ORIGIN SPECIES	Human
ACCESSIONS	BC104432
VECTOR	pDONR223
NCBI SYMBOL	HCG27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagacataccaggctgggcagctgtgggggtccttcagagaggagctgagatacgctacctggaggggccccctgggctggaggggctcctcagtgactgggtgaagtgtttcagaggaccagggtgaggttgggggcatctcatccagaccctgccggcatctgcccc
agaaccaaggggccccctcctcctcctcctcagtggaatgtggggatgtcctcagtcaccccttgagcactcacatcaccccttatgtgaaattttctcactctaa