

C12orf67-chromosome 12 open reading frame 67 Gene

Reference: PTXBC104436

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C12orf67; FAM191B; long intergenic non-protein coding RNA 477		
For research use only.			

GENE DESCRIPTION	chromosome 12 open reading frame 67
GENE ID	144360
ORIGIN SPECIES	Human
ACCESSIONS	BC104436
VECTOR	pDONR223
NCBI SYMBOL	C12orf67
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgccctttttccaacaccacccaagtcggtgtctgtttctcctttcaaggggtctctgctaccccaactctctttctttttttcttctcttgcagagccggctcctcaatgactggctgctttaccttcttcttggaactcattttttcttcgccccgttttagggccatctctatggggatg
tatctcggagcctccactttgactggcttcttttactcgtttcttagggcagctcctcatggacctggagggttggagtgcttggaggggctcctcctctggtggattttttttcttcttcgcttcccaccgtgttccctggggctctgtgtctcgtcttgcctcaaggctacctgttctc
accctcttccccactgagcaagtgctccacagcatccccgacctcgcagaaccaagggtcctggactgagaaagctgtgggtctctgggccccctttttcttcttaa