

## C11orf72-chromosome 11 open reading frame 72 Gene

Reference: PTXBC104448

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C-K-RAS; CFC2; K-RAS2A; K-RAS2B; K-RAS4A; K-RAS4B; KI-RAS; KRAS1; KRAS2; NS3; RALD; RASK2; c-Ki-ras2; GTPase KRas; K-Ras 2; K-ras p21 protein; Kirsten rat sarcoma viral oncogene homolog; PR310 c-K-ras oncogene; c-Ki-ras; c-Kirsten-ras protein; cellular c-		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	chromosome 11 open reading frame 72
<b>GENE ID</b>	283135
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC104448
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	C11orf72
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon )

atgacacaactgccagaactcggcctgaggtctctcaacaacaatccccactggtccccaccactggagcaactctgccaggcttcaagcgaagaagagatccactctgatgcatctcccaggaccttctctgcageatctcaggtcctgctcaccctcagaccat  
cccctctgtgtaactctgtatatacagctcccagccagctccagccacaagccaagaggtatcaaggctctccctgcccactggccccacatccaactgtgtctatctggtcaactttacctccaacattccctgcacagtgatagatggctgtggccctgtcatctcactgaaa  
ttaacaatgatccttggggcctcccaccctacatctggctctcagccccctctctgcaaacatggaacagtgggattatataaatcagaccgttttggcccccttctgctgagagtttttgggtlccccactgccctcagaacaagctcaagacccttggcctggcattcaagtc  
ctttgtctgccactgtcaacctcccagttccagctccttggcttggccatacaagctggtatgcaccagagatctgtctctccaacttggccctagtctctaggagcccttggagctcacctccccaggagcttctgcaaccacctacaacctgtccctataacttataa  
aaccaattag