

C16orf54-chromosome 16 open reading frame 54 Gene

Reference: PTXBC104466

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein C16orf54; chromosome 16 open reading frame 54		
For research use only.			

GENE DESCRIPTION	chromosome 16 open reading frame 54
GENE ID	283897
ORIGIN SPECIES	Human
ACCESSIONS	BC104466
VECTOR	pDONR223
NCBI SYMBOL	C16orf54
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgttgactccagagccgcctctggcgcgctggaggggccccgcacatgggaagcagccccatggccctcactgccctgtgggcctgcaccccatcatgctgtgctggccaccctggctgcgtcttcatcctaccacgcgtgtgttgctgaacgcctgttcgccgtgctc
tcgccagagccccagccacgtgcacccacctgggtgtggcgccagagggagagctgtggaattgagccatgggcacccgcccagagcgtctgaggactggatggctctgcgggtcccccctgctgacagatcggccccctgagcctcccacccaggtgggactttggagggc
cgagcaacagccccacgtgccccctcagccccaaattctctccagcaacttggccccagaccgtactggaggtccagcccgagcaccttctggggccccagccctggaggggagggccccccgcccagggcctggtgagctgggctgaacccgagcagaggccagag
gccagcgtccagtttgggagccccagggcagggcagcggcagggagcccgatcctgagtggggcctccagccacgggtcaccttggagcagatctcagcttctggaagcgtgaaggccggaccagtgtgggttctga