

ODF3L2-outer dense fiber of sperm tails 3-like 2 Gene

Reference: PTXBC104469

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf19; outer dense fiber protein 3-like protein 2; outer dense fiber of sperm tails 3 like 2		

For research use only.

GENE DESCRIPTION	outer dense fiber of sperm tails 3-like 2
GENE ID	284451
ORIGIN SPECIES	Human
ACCESSIONS	BC104469
VECTOR	pDONR223
NCBI SYMBOL	ODF3L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggggaccctcagctgcactccacccacggctggccacagcccccttggccgagtgacggagggccagattccggagacggcctgaggaagtctgtgggacggccacccctggagaacggctggggccaggttgcctgctgcctccacgtgggtcatcaaccacgactgcaccaggggtggccagtcccgctactcgtcgtccggagggccagtgaggcaccacctcaggacaccagccggggcccatctacttcttgaccccaagtcacacgttggccgagctgcacccctgcctactccatgcaggccggccaaagtctcgggtccggaggtgacgccaggccctggggcctacagcccagagaaggtgccccctgcgcgccatcggacaccgcccgtttcacctgggctgcgcctccccctgaagccctggacacctcagccccctgccccctaatgcctacaccatgccccctctgtggggctcacagatcttaccgaagcccagcagccccagctacacgggtggggccgaccccccccgcccccccgagatctgtctgagataccagggccggccagtagcacagcccgagcgaacacacctaccgcccagcgctgctgccttaccatgtggggcgccccggccccggcaccctggaaggacccccggccctgcgcccactgcccagagcaggtcacctggaacaaagccagggtccagcgttctctatggccatccgcactcgaacggggccagcaccatggccgccaccagccctcccgccctcgggggcacaggctgcccggccgtgctgtag