

LOC100049076-similar to SMA4 Gene

Reference: PTXBC104473

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to SMA4
GENE ID	100049076
ORIGIN SPECIES	Human
ACCESSIONS	BC104473
VECTOR	pDONR223
NCBI SYMBOL	LOC100049076
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgattgctcaccaaaagcctggaccctcccagcctgtgaccttgggaccaactccactacgcagcagacaaggggctctgtatgtggatgtgatccgtgtgaacagctactacttggatcgcactacgggcacctggagttgattcggtgcagctggccgccagttt
gagaattggtgaagacatcacatccattatcagagcgcgtatggagtggaaacgctttagggcttcaccaggttaagcgtgttga