

MCM3AP-minichromosome maintenance complex component 3 associated protein Gene

Reference: PTXBC104958

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAP80; SAC3; germinal-center associated nuclear protein; 80 kDa MCM3-associated protein; MCM3 acetylating protein; MCM3 acetyltransferase; MCM3 import protein; MCM3 minichromosome maintenance deficient 3 associated protein; germinal center-associated nucl		

For research use only.

GENE DESCRIPTION	minichromosome maintenance complex component 3 associated protein
GENE ID	8888
ORIGIN SPECIES	Human
ACCESSIONS	BC104958
VECTOR	pDONR223
NCBI SYMBOL	MCM3AP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr