



Reference: PTXBC104992

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1033; MRT43; SWIP; WASH complex subunit SWIP; WASH complex subunit 7; strumpellin and WASH-interacting protein; WASH complex subunit 4		

For research use only.

GENE DESCRIPTION	KIAA1033
GENE ID	23325
ORIGIN SPECIES	Human
ACCESSIONS	BC104992
VECTOR	pDONR223
NCBI SYMBOL	KIAA1033
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]