

ARL16-ADP-ribosylation factor-like 16 Gene

Reference: PTXBC105076

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ADP-ribosylation factor-like protein 16; ADP ribosylation factor like GTPase 16		
For research use only.			

GENE DESCRIPTION	ADP-ribosylation factor-like 16
GENE ID	339231
ORIGIN SPECIES	Human
ACCESSIONS	BC105076
VECTOR	pDONR223
NCBI SYMBOL	ARL16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagagtggcgggtggcgccgtgagccggcgccgagctacgggtgcccgggtggagcgaagcacggaatgtgtctctgctggggccacggcgctgggaagacgctgctgggaaacggctgcaggaggtgagctccccgggatgggaaagcgacgtgggggagc
cgccccgacacggcccggtggcgacacatcttactgacatctgtggcacagagaagatcaccatccgggagcttgggggtgcatggggcccatctgtccagttactatggaactgccgttctctctgtttgtgatggagcctctgacccaccagctctctgcatctgtgtg
cagctcttaggtctctcttctgcagaacaactgcagaagcctcgtgctgatactcttcaataaaatcgacctaccctgttacctgtccacggagagatgaagtcattaatcaggcttcagacatcattgcttggccaagcagaacatcaccacggcagaaatcagcggccgtgaaggca
ctggcttagcagggtgctggcctggctccaggccaccacagagccaacgattga