

LOC339788-hypothetical LOC339788 Gene

Reference: PTXBC105688

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized LOC339788		

For research use only.

GENE DESCRIPTION	hypothetical LOC339788
GENE ID	339788
ORIGIN SPECIES	Human
ACCESSIONS	BC105688
VECTOR	pDONR223
NCBI SYMBOL	LOC339788
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcccaggagctgcctatggattcaagaaaactgtatgagaggagtaacatttgggtggaagcgaatacaagtaaacatatgtcaaatgggtctggaagcaatcctccacagaaagaaagcagaatcttctgggcatttctgagcaagtggactcctcagctgagcagctc
ctcaggagaagcggacagggtgcaagtggtggaagccaacactcacagcagttgggtggtgagtgctctggcttccaacagcattgaaaccaacagcattgactacagcgagttctcaatctcagaacacatcccagcctccctggggctgctaa