

LOC203547-hypothetical protein LOC203547 Gene

Reference: PTXBC105694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein LOC203547
GENE ID	203547
ORIGIN SPECIES	Human
ACCESSIONS	BC105694
VECTOR	pDONR223
NCBI SYMBOL	LOC203547
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagcggcgataaggcggcgctgaacgcactgcagcctcctgagttcagaaatgaaagctcattagcatctacactgaagacgctcctgtttctcacagctttaatgatcactgttctattgggttatattcacactaaatcttacatattgaaggcgcccttgggatgtccaatagg
acagctattttacgctgctattgttcagtggtcgcctccatgtgtgtccttctgtgtatgtggcctggaatgaaggctcagacagtggcgtgaaggcaaacaggattaa