

LOC91948-hypothetical LOC91948 Gene

Reference: PTXBC105716

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical LOC91948
GENE ID	91948
ORIGIN SPECIES	Human
ACCESSIONS	BC105716
VECTOR	pDONR223
NCBI SYMBOL	LOC91948
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgcggcctgccaggagagaataaacagcaccggcagcctggcacagagctctggcgggaacgtccaccaagatcggtcaggcactctcatggcgtcctcctggatacaggcttttgcctttgttttaatttggagcagtgaaagccgaaagtgaggagtgaggacaaggagaag
accatttgcctctgttcaatcaaacagcatctctcatttgaattggcaagaggcagttcggggctgtggcaggctctgtacaccacctgtgcctccaagtcactgtgctttctggtttctttacctatgtctgtaaaacgccagcaagaactctacatctcaacattttag