

## ZNF816A-zinc finger protein 816A Gene

Reference: PTXBC105930

|                               |                                                            |                     |                |
|-------------------------------|------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                        | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                       | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | ZNF816A; zinc finger protein 816; zinc finger protein 816A |                     |                |
| <b>For research use only.</b> |                                                            |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | zinc finger protein 816A                                                                         |
| <b>GENE ID</b>          | 125893                                                                                           |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC105930                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | ZNF816A                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgttacgtgaggagccaccaagaagcaaaaggagccaggatgctctctcaggagacgttgaacttcaggatgtgctatagagttctcttggaggagtgaaatgcctgaacctgcacagaggcttatacagggtgtgatgttgagaactacaggaacctgg  
agtttggatagctctttaaataccatgatggagttctcatcaaccaggccacagtattacaggagaagtgtatccacagggacgttgcaagacataaaagtcacacattggagatttttctccagaaatgaagaaatattcatcatttgagttcagtgccagaaagtgaagaa  
atggccatgaagcaccatgacaaaaatcaaaaagtgtgactgttagtacagaccgaagtgtatcacagggcatgctggaacaagcctattaaagatcagcttgattagcttctcgcacatcgcctgaactccacatgtttcagactaaaggaaaattgtaaccaattggacaagtctatc  
ggctctcctcagcttcagaatcccaagaatttctgttaggtcacaactcatatttctaataagatgggaagaatttctccattcttcattcacacaatacaggaataatgcatgagagaaaaaccttgcacaagtaagtgtgtggcaagccttaattatagctcactcttaaggagacac  
cacataaccatcaagagagagagataataatgtgatgtgtgcaagatctttaatcagaagcaatacattgtatatcatcatagatgtcacactggtgagaaaaactacaagtgtatgtgtgggaagaccttcactcagatgtcatccttgtatgccatctgtagactcactggag  
agaaaccttacaagtgtatgtgtgtggcaagaccttcagtgagaaatcatccttagatgccatctgtagactcactcagagagaaaccttacaagtgtatgtgtgtggcaagaccttggcgaattcagccttgaattcataaggcaattcactcagagagaaaccttacaagtgtat  
atgagtggtggcaagaccttcagtcagaatcatccttcaatgccatcatatcttcacactggagagaaaccttacaagtgtgagaatgtgacaatgtttacattcgagatcacaccttgaaagacatagggaagattcactcagagaggatcatcaaatgttaaggtttgtgacaaggtttt  
ccggagtgattcactcctgcagaacatcagagagttcactcggagaaaccatacaagtgtataaatgtggcaggagtttcagtcggaagtcatccttcaataccatcacacttcacactggagagaaaccttacacatgtaagtgtggcaaggtttttatgtcgaagagaaaaac  
ttgcacgtcatagactcagctgtgagagaaaccttacaagtgtgagaatgtgacaagtgtttcagtcgagatcacaccttgaagacataggagaattcactcgggaaaaaccatacaaatgtaaggtttgtgacaaggcttccggagtgattcagccttgcacaaccatcagaga  
gttcatactggagagaaaccttacaagtgtataaatgtgcaaggtttttatcaaaaaggaaatccttcacaacatcagagagttcactcggagagaaaccttacaagtgtatgaaatgtggaaggtttttatcaaaaaggcaagccttgcacaacatcagagagttcactcagagaa  
accttacaagtgtatgtgtggcaagccttactggacagtcacacttattcaccatcaagcaatccatgggtgtagggaaccttacaagttaa