

NEURL2-neuralized homolog 2 (Drosophila) Gene

Reference: PTXBC105935

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf163; OZZ; OZZ-E3; neuralized-like protein 2; neuralized homolog 2; neuralized E3 ubiquitin protein ligase 2		

For research use only.

GENE DESCRIPTION	neuralized homolog 2 (Drosophila)
GENE ID	140825
ORIGIN SPECIES	Human
ACCESSIONS	BC105935
VECTOR	pDONR223
NCBI SYMBOL	NEURL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtgctcgcagcccggtggtgactctgggactcagcgcggagcccccctccaccgcttcacgtgacgtgccaacatccggtggaccctctgggacgcgggcccacacgcgtggagagcttgcgccacggcgtgtgttcagcgcgcgagcc
gctggcccccggccagggtcttctgtcgaatcgaggagaaagagctgggtgtgtgacacatctgctcgtgtgacgcgtggaccccgccagctctgccccgttcccgagttttctctgcccgatctgttaacctgggcccacacctgggtcttcgccatcagcgccacca
caaccgctgccccgggagggccgcccggagggcagcggccccacgcgcacccaacctcctctgtggaaccatactgcgcattgagcagtttgcattccccgggacgcctgggtggccgcagccggccagggctctacagccatctctggaccagctctatgagct
gaacgtgctgcctcgcacgcgcgcgttagccgcctgtgttcttttccccgcggccgatggcagggccgacatgcacatcatcaacggcgaggacatggcccgagcggccgggactgcccagctgcgcagccctctacgggtgtgtggacgtgtttgcttcacaaa
gagcgtgctgcttgcagctcagatggtgtgcatccctgcagactctgtccgcttagtgatacaaaaggagcatggtgcaccgctggtgacattgaggtccacctgcccagaagaacttaaggatttctgcaagtatgagtgatga