

ZDHHC8P-zinc finger, DHHC-type containing 8 pseudogene Gene

Reference: PTXBC105936

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	zinc finger, DHHC-type containing 8 pseudogene
GENE ID	150244
ORIGIN SPECIES	Human
ACCESSIONS	BC105936
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC8P
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatccatccaggagcgaaggacaggaggagcgtgagtgctgctgcctccaggccgactcactcttcggcgactcaggcgtctatgatgctcccagctctacagcctgcagcaggccagtgctgctgagggcttcgaggtcctacgtgtgctacagctctacagatgaccttgtggccaggcccgcttcggcgccgctgcaaccctgtcctgcagacatcattgtcctcgtttccagctccgtgagccgtgcactgaggacgtcgtcctccctgcaggctgacaggtggggaagctgaggcaggaagccctttag