

BTG2-BTG family, member 2 Gene

Reference: PTXBC105949

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein BTG2; APRO1; PC3; TIS21; B-cell translocation gene 2; BTG family member 2; NGF-inducible anti-proliferative protein PC3; nerve growth factor-inducible anti-proliferative; pheochromocytoma cell-3; BTG anti-proliferation factor 2		

For research use only.

GENE DESCRIPTION	BTG family, member 2
GENE ID	7832
ORIGIN SPECIES	Human
ACCESSIONS	BC105949
VECTOR	pDONR223
NCBI SYMBOL	BTG2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccacgggaagggaaccgacatgctcccgagatgccgcccgtgggtctctccagcctctgaggacccgggctgctgagcagcagaggettaaggtcttcagcggggcgtccaggaggcactcacagagcactacaacaccactgggtccgaaaagccgtccaagggtccgctaccgctgcatcgcatcaaccacaagatggaccccatcatcagcagggtggccagccagatcggaactcagcagcccgagctgcaccagctgtgccagcagctgacccctgtgggtggaccctatgaggtgtcctaccgcatggggagggagggctcctctgctgtgtacgaggaggccctgcccgcctctgtgggtctctacctgcaagaaccaagtgctgctgggcccggagcagccctccaagaactacgtgatggcagctctcagctag