



Reference: PTXBC105960

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SP140 nuclear body protein; nuclear body protein SP140; LYSP100; LYSP100-A; LYSP100-B; lymphoid-restricted homolog of Sp100; lymphoid-specific SP100 homolog; nuclear autoantigen Sp-140; speckled 140 kDa		

For research use only.

GENE DESCRIPTION	SP140 nuclear body protein
GENE ID	11262
ORIGIN SPECIES	Human
ACCESSIONS	BC105960
VECTOR	pDONR223
NCBI SYMBOL	SP140
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr