

INSL3-insulin-like 3 (Leydig cell) Gene

Reference: PTXBC106722

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	prepro-INSL3; RLF; RLNL; ley-I-L; insulin-like 3; insulin-like 3 (Leydig cell); leydig insulin -like hormone; leydig insulin-like peptide; relaxin-like factor b; insulin like 3		

For research use only.

GENE DESCRIPTION	insulin-like 3 (Leydig cell)
GENE ID	3640
ORIGIN SPECIES	Human
ACCESSIONS	BC106722
VECTOR	pDONR223
NCBI SYMBOL	INSL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacccccgtctccccgctggcgctggtgctgtggccctgccctggtgttcggtggccccgcgccccagagatgcgtgagaagttgtgcggccaccacttcgtacgcgcgtagtcgcgctgtgcggggggccccgcgtgtccaccgaagccaggaggcctgc
gaccggaggcgaccgtgagttgctacagtggtggagagacgacatctgctccatgggctggtggccgacagtaatctcacgtgggacctggcctgcagccctgccccagaccttcacatcacgcccaccacgtgcagctgccaccaacctgcacgctactgctgcctcagt
ggctgtaccaacaagacctgctgacctctgtccctactga