



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CEACAM3-carcinoembryonic antigen-related cell adhesion molecule 3 Gene

Reference: PTXBC106728

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CD66D; CEA; CGM1; W264; W282; carcinoembryonic antigen-related cell adhesion molecule 3; CD66d antigen; carcinoembryonic antigen CGM1; carcinoembryonic antigen gene family member 1; nonspecific cross-reacting antigen; carcinoembryonic antigen related cell

For research use only.

GENE DESCRIPTION	carcinoembryonic antigen-related cell adhesion molecule 3
GENE ID	1084
ORIGIN SPECIES	Human
ACCESSIONS	BC106728
VECTOR	pDONR223
NCBI SYMBOL	CEACAM3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggccccctcagcctctccccacagagaatgcatccctggcaggggcttctgctcagagcctcacttctaactctfaggaacggccaccactgccaaagctcactatgaaatcatgcccgtcagctgcagaggggaagagtgctctcactgtcccaactctccccagca
tcttttggctcagctgggtacaaaggggaagagtggtggcacaagctctaattgtaggatgtaataaggaaactcaaacagctaccccaggggccgcatacagcgctgcagagacaatataccaaatgcatccctgctgatcagaatgtcaccagaaatgacataggaattctacacc
tacaagtcataaagtcagatctgttgaaagaaagcaactggcagcttcatcgtatatacagaanaatgccccaggccctctctgtggggcgctgcgcgcgcagctgaccgggggtctgctggagtgggcctgtgtggccgcgtgtgtgtttctctctcttgcacaaactgcgaagaac
cagatccagcgtgactctcaaggagcagcagccccaaagcccttgccttggcctgtgctccctggatcagctcctcctcactgccagggcccccctaccacacccaggacagcagcttccatctatgagaattgctaaacatgacacaaacattactgcc

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1