

GABARAP-GABA(A) receptor-associated protein Gene

Reference: PTXBC106749

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GABARAP-a; ATG8A; MM46; gamma-aminobutyric acid receptor-associated protein; GABA(A) receptor-associated protein; GABA type A receptor-associated protein		

For research use only.

GENE DESCRIPTION	GABA(A) receptor-associated protein
GENE ID	11337
ORIGIN SPECIES	Human
ACCESSIONS	BC106749
VECTOR	pDONR223
NCBI SYMBOL	GABARAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagttcgtgtacaaagaagacatccgttcgagaagcgcgctctgagggcgagaagatccgaagaatacccgaccgggtgccgtgatagtagaaaggctccaaagctcgataggagacctggacaaaaagaataacgtggtgcctctgatctcacagttggtcagtt
ctacttctgatccggaagcgaatcatctccgagctgaggtgctgtttttctttgtcaacaatgtcatccaccacagtgccacaatgggtcagctgtaccaggaacacatgaagaagactctttctacattgcctacagtgacgaagtgctctacgggtctgtga