



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SLMO1-slowmo homolog 1 (Drosophila) Gene

Reference: PTXBC106750

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SLMO1; C18orf43; HFL-EDDG1; PRELI domain containing protein 3A; erythroid differentiation and denucleation factor 1; protein slowmo homolog 1; slowmo homolog 1; PRELI domain containing 3A		

For research use only.

GENE DESCRIPTION	slowmo homolog 1 (Drosophila)
GENE ID	10650
ORIGIN SPECIES	Human
ACCESSIONS	BC106750
VECTOR	pDONR223
NCBI SYMBOL	SLMO1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatctgaggcctcgagcacgctgtttggccaccctggggacacggctcatccagcggccatgcgcgaagtaccgcaacccgatgaacccgagcgtcgtcggcgctgagtgctctacagcgccgcgtggagccgcgcgcctgcacagcttgcgcctgctcagaccagtg
ggggctggccagctctcgtgagcagcatttggccacaggtaggacatgcatcatcaggaacatcttcggatggatccgtggaaagaaatggcaacttggctcaacaatfacacacacacaaatttgcgtgattagagaggtgtgtgtaacacatctcatcaggaacccagaa
atacctctctcagcaagaagccatcatcactctctcctgctgctgattcgtttgaaatttgaatggcctcaatcatalcctaactcaaaaagggctggcctctatcagtcgaatttgaacatctctaaacgcctctctgaacta

More info: Please contact us for any question: contact@proteogenix.fr