



Reference: PTXBC106752

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD37 molecule; leukocyte surface antigen CD37; leukocyte antigen CD37; GP52-40; TSPAN26; cell differentiation antigen 37; tetraspanin-26; tspan-26		

For research use only.

GENE DESCRIPTION	CD37 molecule
GENE ID	951
ORIGIN SPECIES	Human
ACCESSIONS	BC106752
VECTOR	pDONR223
NCBI SYMBOL	CD37
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]