

FBXL21-F-box and leucine-rich repeat protein 21 Gene

Reference: PTXBC106753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBL3B; FBXL3B; FBXL3P; Fbl21; F-box/LRR-repeat protein 21; F-box and leucine-rich repeat protein 21; F-box and leucine-rich repeat protein 3 pseudogene; F-box and leucine-rich repeat protein 3B; F-box protein Fbl3b; F-box/LRR-repeat protein 3B; F-box and		

For research use only.

GENE DESCRIPTION	F-box and leucine-rich repeat protein 21
GENE ID	26223
ORIGIN SPECIES	Human
ACCESSIONS	BC106753
VECTOR	pDONR223
NCBI SYMBOL	FBXL21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcactcttctatatgaagaggaattcgagacgttctcaagaagaacccctgttactcactttatttggctcagtcagcaaatgggttttaggacgggtaggctcaactgtcctcgactgattgagttagtggtgtgtgctaatgatcttcagcctcttgataatgaactatttattgctgaacactgtacaaacctaacagccttggccctcagcaaatgtgaagttagctgcagtcctcatcaggtttgaagactgtgtgagagaagggttaacacagctctctgtaattggaggaaatttggatccctgatgaggattatagcctagatgaaattcacactgaagctcctcaatacctggaagagtatggctcctgatgtgatgcctctctggttaa