

HCP5-HLA complex P5 Gene

Reference: PTXBC106759

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	6S2650E; D6S2650E; P5-1; HLA complex P5 (non-protein coding)		
For research use only.			

GENE DESCRIPTION	HLA complex P5
GENE ID	10866
ORIGIN SPECIES	Human
ACCESSIONS	BC106759
VECTOR	pDONR223
NCBI SYMBOL	HCP5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctgctgaggatgtcagagcacaggagccttgggaattacctggaaatgcgactgaaatctctctgaggggtctgggtcttggaaatcaaaccctctcaggttgggtggctggacgattctctcacactacaatgggacaaggggaaccaggaggccccaaggg
gatccctgggttccacacgaactctctacacctcattgtgtgacagcagccatgcctctctctggggatcaggatctattacctgtgcctggagaggaggggactctctctcaccctgggtctctggacacatactgtccaattccctgtggcagctgtaatgtgtgttaatgggcact
cattgtcccttttaagggtaccctccttag