

## ETV1-ets variant 1 Gene

Reference: PTXBC106762

|                               |                                                                                        |                     |                |
|-------------------------------|----------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | ETS translocation variant 1; ets variant gene 1; ets-related protein 81; ETS variant 1 |                     |                |
| <b>For research use only.</b> |                                                                                        |                     |                |

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ets variant 1                                                                                    |
| <b>GENE ID</b>          | 2115                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC106762                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | ETV1                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atggatgatttaccagcaagtgccttacatggtcaccaatagtcagcgtgggagaattgtaacgagaacacaaatgtcaggaaaagaaattcattaacagagatctggctcatgattcagaagaactcttcaagatcaagtcaattacaggaaacatggcttcagaagctc  
aggtacctgacaatgatgagcagtttaccagactatcaggctgaagtgttgccttcatggcctgcaactgaaaatcaagaagaacccacagtcacatgttcagaatcagctcgtcagtcagaacagcccttaattcagctatggagaaaagtgcctgtacaatgtcagtgcc  
tatgatcagaagccacaagtgggaatgagccctcaacccccacaccatcagcagcaggtgtccctcagtcacatgcatctccaaactcaactcacaacgaaactgaccggcctccagctcactccctccatcagtcacatcagatagcagctacccatggac  
cacagatttcgccagccttctgaacctgttaactccttctccttgcgacgatgccaagggaaggacgtcctatgtaccaacgccagatgtctgagccaacatccctccaccacaaggctttaagcaggagtagccagaccagctgtatgaaccaaacacatggttggcagt  
gcggccagccaaagcttccctcctctgatgattaaacaggaacccagagatttgcatactcagaagtgctgacgtccactccatttatagggcaagaaggcttcctgctcatccagcagaacagaaggctgtatgttgaaggcccccaggcagttttatgatgacacct  
gtgttgcagaaaaatcagtgagacatcaaaagagccaggaatgtatcgggaaggaccacataccaacggcgagatcactcagctctggcagtttttgtagctcttctgatgaccttcaaatctcattttatgcctggactggtcagggatggaatttaactgattgagc  
ctgaagaggtggccgacgttggcattcagaaaacagccagctatgaactatgataaacttagccgttactccgtattactatgagaaggaaattatgcaaaagggtgctggagagagatgtctacaagtttgggtgatccagaagccctttctcatggcctttccagataat  
cagcgtccactgctgaagacagacatggaacgtcacatcaacgaggaggacacagtcctcttctcatttgcagagagatggcctacatgccggaagggtgctgcaacccccaccctacaacgaaggctacgtgtattaa