

## TCTEX1D1-Tctex1 domain containing 1 Gene

Reference: PTXBC106885

|                               |                                                                |                     |                |
|-------------------------------|----------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                            | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | tctex1 domain-containing protein 1; Tctex1 domain containing 1 |                     |                |
| <b>For research use only.</b> |                                                                |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | Tctex1 domain containing 1                                                                        |
| <b>GENE ID</b>          | 200132                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC106885                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | TCTEX1D1                                                                                          |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgatgatgcagacaatgctaaaggcagagcagctcatcatggaagaaaaggaggatattctcttaagtaaatcatgaatttggcgaaggaaatcatggcgcatcaaaagattctatgagtactgtgtcttatatggaagaaccagtcagcgtgatgatctctgccttacagttc  
agatggaaaacacctatcagttgggtctcccaaacatttctgtgtgcaccgtcaatcatatttgaagatgtagtaaccagctatctacaagtagaagaatatgaaccagagctctgtagacagatgactaaaaccatttctgaggttattaaagcccagggtcaaggacttgatgattccacg  
gtataaactaattgtgattgttcacattggacaactgaacaggcagagcataattattggaagcagatgcctctgggacctaagtgatacctttcatcttatgttttcagaaattctctctcttcgctcttgcaaatgctatgcagtttaccttgagtgga