

CPLX4-complexin 4 Gene

Reference: PTXBC106887

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPX-IV; CPXIV; complexin-4; CPX IV; complexin IV; complexin 4		
For research use only.			

GENE DESCRIPTION	complexin 4
GENE ID	339302
ORIGIN SPECIES	Human
ACCESSIONS	BC106887
VECTOR	pDONR223
NCBI SYMBOL	CPLX4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcttcttatgaaaagtatgataagtaaccaggtaaagaatttaggtttggtgggtctgaagaaaataaagaaggaggtgcacatctgcagcagctcaagggtgactagagaggatgaggagatcaaaagcaaatgattgaggagaagatggaaagagatg
ctgcatttacagaaaaaggcagaaaggcatgctcagagttcatctcagagaaaaatacaggctccaaagagtgaaatggatgagaatcaaatccagatggctggagatgattggattacctgaagatctccgaaaaatggtagatgaagatcaagaagggaagataaag
attctattctgggcagatacagaatctccagaactggacttgataccataaaaaaaaagccaggccaccttcaatgaaatcaagcagacagcggagcagaagtgttcctggtgatgga