

SYNGR4-synaptogyrin 4 Gene

Reference: PTXBC106891

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	synaptogyrin-4; testis tissue sperm-binding protein Li 72n; synaptogyrin 4		
For research use only.			

GENE DESCRIPTION	synaptogyrin 4
GENE ID	23546
ORIGIN SPECIES	Human
ACCESSIONS	BC106891
VECTOR	pDONR223
NCBI SYMBOL	SYNGR4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacatccccaaagcctccaggagctggccaacagcgaagccgtgcagttctgagaaggccaagaccatcacgcgggtcttgaagggtcttctccctgatcgtcttctctctctgctgaccgacggctaccagaacaagatggagttcgcgcagctccactgcattctcaaca
gcaacacgctggcctgcagcttgcctggggagcggccttctggccttctcagctgcctggccttctctgctgacacacagagaccgcatggcggcaccgcttcaagacagccttccagctcctggacttcatctggtgttctctggcagttgtctgttcatgggttctg
cttctggccaacaatggcagcattcggcggccaagagttctctctggggagcagcagtgccaggcagccatcgccttcttctccatcttctgtggaattccaggcctacctggcattccaggacctccgaaatgatgctccagctcccttacaagcgttctggtgaggggt
ggcatgtgtgacccaccctcccttgcctctgccaacagccctgtgaacatgccaccactggcccaacagcctgagttatgtagctctgcccctgtctgaccgctccaaagtcccccggttctgtatgatgctgacaactaa