

LRRC10-leucine rich repeat containing 10 Gene

Reference: PTXBC106893

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LRRC10A; leucine-rich repeat-containing protein 10; leucine rich repeat containing 10		
For research use only.			

GENE DESCRIPTION	leucine rich repeat containing 10
GENE ID	376132
ORIGIN SPECIES	Human
ACCESSIONS	BC106893
VECTOR	pDONR223
NCBI SYMBOL	LRRC10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgtggacaagatggatctgagtgaggagccagttacgcggttccccctgcacgtgtgtctctcagggagctgtcaagctctacctgagcgacaaccacctaatagcctgcctcggagctggggcagctacagaacctgcagatttggccttgatttcaacaactcaa
ggctctgccccagggtgtgtgacattgaacagctctgcatctctacctgggcaacaacaactctgcgacctccccagtgagctgagcctgtccagaacctcaggacctgtggatcgaggccaactgcctcaccagctgccggatgtgtgtgtgagctgagctccttaagactc
tgcatgccggctccaacgccctgctgttgcagggccagctccggcgctccaggaagctgaggaccatctggctctcgggcaaccggcctaactgacttcccactgtgctgtctcactgccccttctggaggtgattgatgtggactggaacagcatccttactccccagctggcg
cacctgtcaagctgaagctggtcatctatgaccacaatccttgcaggaacgcaccaaggtggccaaggtgtgcgctgtgtgggagatgggcagaggagacgcagagccccgacctagaaaagccaggcgctatgcgttggtcagagaggaaagccaggagctacaggcac
cagtcctctacttctctaccaactcctga